

Code No.13-0306

## SPECIFICATION

## Methyl Salicylate

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ slightly yellow, clear liquid
Solubility in potassium hydroxide solution	to pass test
Density (20°C)	1.180~1.190g/mL
Acidity (as HOC <sub>6</sub> H <sub>4</sub> COOH)	max.0.1%
Heavy metals (as Pb)	max.0.001%
Assay (capillary GC)	min.98.0%