

## SPECIFICATION

## 2-Hexanone

Wako 1st Grade

| REQUIREMENT                       | SPECIFICATION                             |
|-----------------------------------|---|
| Appearance                        | Colorless ~ slightly yellow, clear liquid |
| Density (20°C)                    | 0.807~0.814g/mL                           |
| Refractive index 20°C             | 1.398~1.402                               |
| Water                             | max.0.8%                                  |
| Acidity (as CH <sub>3</sub> COOH) | max.0.05%                                 |
| Assay (capillary GC)              | min.95.0%                                 |