

Code No.13-0406

SPECIFICATION

Methyl Decanoate

Wako 1st Grade

| REQUIREMENT | SPECIFICATION |
|-------------------------------------------------------------------|--------------------------------------------|
| Appearance | Colorless ~ nearly colorless, clear liquid |
| Density (20°C) | 0.869~0.876g/ml |
| Refractive index 20°C | 1.424~1.429 |
| Water | max.0.2% |
| Acidity [as CH ₃ (CH ₂) ₈ COOH] | max.0.1% |
| Assay (capillary GC) | min.95.0% |