

SPECIFICATION

4-Methyl-2-pentanone

for Atomic Absorption Spectrochemical Analysis

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Water	max.0.05%
Copper (Cu)	max.2ppb
Zinc (Zn)	max.5ppb
Cadmium (Cd)	max.2ppb
Lead (Pb)	max.1ppb
Chromium (Cr)	max.5ppb
Manganese (Mn)	max.1ppb
Iron (Fe)	max.5ppb
Cobalt (Co)	max.2ppb
Nickel (Ni)	max.5ppb
Peroxides	to pass test
Assay (capillary GC)	min.99.5%