

Code No.13-0776

## SPECIFICATION

## Methoxyacetic Acid

Wako 1st Grade

| REQUIREMENT                 | SPECIFICATION                             |
|-----------------------------|---|
| Appearance                  | Colorless ~ slightly yellow, clear liquid |
| Solubility in water         | to pass test                              |
| Density (20°C)              | 1.172~1.180g/ml                           |
| Refractive index 20°C       | 1.412~1.418                               |
| Water                       | max.3.0%                                  |
| Assay (volumetric analysis) | min.95.0%                                 |
| Assay (GC)                  | min.95.0%                                 |