

Code No.13-0849

SPECIFICATION

trans-2-Methyl-2-butenal

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ greenish yellow, clear liquid
Solubility in ethanol	to pass test
Density (20°C)	0.864~0.871g/ml
Refractive index 20°C	1.443~1.449
Water	max.1.0%
Acidity [as CH ₃ CH:C(CH ₃)COOH]	max.2.0%
Assay (capillary GC)	min.97.0%