

Code No.13-0882

# SPECIFICATION

## Methanol 300

for Pesticide Residue-PCB Analysis

| REQUIREMENT                                    | SPECIFICATION          |
|--|------------------------|
| Appearance                                     | Colorless clear liquid |
| Density (20°C)                                 | 0.790~0.793g/ml        |
| Water  | max.0.05%              |
| Residue after evaporation                      | max.1ppm               |
| Suitability for pesticide residue-PCB analysis | to pass test           |
| Assay (capillary GC)                           | min.99.8%              |