

## SPECIFICATION

## Methanol

for Preparative Chromatography

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Density (20°C)	0.790~0.793g/ml
Refractive index 20°C	1.327~1.330
Absorbance (220nm)	max.0.30
Absorbance (260nm)	max.0.02
Absorbance (280~400nm)	max.0.01
Water	max.0.05%
Residue after evaporation	max.5ppm
Acidity(as HCOOH)	max.0.001%
Peroxides (as H <sub>2</sub> O <sub>2</sub> )	max.5ppm
Assay (capillary GC)	min.99.7%