

Code No.13-1233

SPECIFICATION

Magnesium Sulfate

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White powder
Solubility in water	to pass test
pH of a 50g/l solution at 25°C	6.0~9.0
Chloride (Cl)	max.0.002%
Nitrate (NO ₃)	max.0.002%
Heavy metals (as Pb)	max.0.001%
Sodium (Na)	max.0.01%
Potassium (K)	max.0.01%
Calcium (Ca)	under investigation
Arsenic (As)	max.1ppm
Iron (Fe)	max.0.001%
Ammonium (NH ₄)	max.0.003%
Assay	min.98.0%