

SPECIFICATION

Metaphosphoric Acid, Lump

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless or white mass
Solubility in water	to pass test
Chloride (Cl)	max.0.002%
Nitrate	to pass test (as NO ₃ :max.about 0.001%)
Sulfate (SO ₄)	max.0.01%
Calcium (Ca)	max.0.01%
Lead (Pb)	max.0.002%
Arsenic (As)	max.3ppm
Iron (Fe)	max.0.005%
Substances reducing permanganate	to pass test (as H ₃ PO ₃ :max.about 0.02%)
Assay (as HPO ₃)	min.37.0%