

SPECIFICATION

Magnesium Nitrate Solution (10mg Mg/mL 15% Nitric Acid Solution)

for Atomic Absorption Spectrochemical Analysis

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Acidity (HNO ₃)	under investigation
Silver (Ag)	max.50ppb
Cadmium (Cd)	max.50ppb
Thallium (Tl)	max.100ppb
Lead (Pb)	max.100ppb
Arsenic (As)	max.100ppb
Selenium (Se)	max.100ppb
Concentration	9.50~10.50mg Mg/mL