

Code No.13-1600

# SPECIFICATION

## Metalaxyl Reference Material

Traceable Reference Material

REQUIREMENT	SPECIFICATION
Appearance	White, crystalline powder ~ powder
Solubility in acetone	to pass test
Melting point	72~75°C
Specific rotation [ $\alpha$ ]D20(c=1,CH3OH) (after drying)	-0.5~+0.5°
Loss on drying over silica gel in vacuum	max.0.3%
Assay	min.98.0%
Assay (capillary GC)	min.98.0%