

SPECIFICATION

Multielement Standard Solution W-XI
for ICP-MS Analysis

REQUIREMENT	SPECIFICATION
Appearance	Slightly red ~ pale brown or pale yellowish brown, clear liquid
Acidity (HNO ₃)	0.45~0.55mol/l
Beryllium (Be)	95.0~105.0mg/l
Cobalt (Co)	95.0~105.0mg/l
Gallium (Ga)	95.0~105.0mg/l
Indium (In)	95.0~105.0mg/l
Thallium (Tl)	95.0~105.0mg/l
Yttrium (Y)	95.0~105.0mg/l