

Code No.13-1644

# SPECIFICATION

## Methyl Benzoate

for the Japanese Pharmacopoeia General Tests

| REQUIREMENT                       | SPECIFICATION          |
|-----------------------------------|------------------------|
| Description (JP method)           | Colorless clear liquid |
| Refractive index 20°C (JP method) | 1.515~1.520            |
| Sp. Gr. d 20/20 (JP method)       | 1.087~1.095            |
| Acid value (JP method)            | max.0.5                |
| Purity (JP method)                | to pass test           |
| Water                             | max.0.5%               |