

## SPECIFICATION

Methanol  
for Medicine Test

REQUIREMENT	SPECIFICATION
* Appearance	Colorless clear liquid
Assay (CH <sub>3</sub> OH) (GC)	min.99.8%(mass/mass)
Solubility in water	to pass test
Density (20°C)	0.791~0.793g/ml
Refractive index 20°C	1.327~1.330
Water	max.0.1%(mass/mass)
Residue after evaporation	max.5ppm(mass/mass)
Acidity (as HCOOH)	max.0.002%(mass/mass)
Alkalinity (as NH <sub>3</sub> )	max.3ppm(mass/mass)
Carbonyl compounds (as CO)	max.0.005%(mass/mass)
Ethanol (GC)	max.0.05%(mass/mass)
Substances darkened by sulfuric acid	to pass test
Substances reducing permanganate (as O)	max.5ppm(mass/mass)
Sp. Gr.(20/20°C) (Ph.Eur.)	0.791~0.793

JIS K 8891:2006(\*:Additional test performed by Wako but not required by JIS)(Ph.Eur.)