

SPECIFICATION

Methanol
for Medicine Test

REQUIREMENT	SPECIFICATION
* Appearance Assay (CH ₃ OH) (GC) Solubility in water Density (20°C) Refractive index 20°C Water Residue after evaporation Acidity (as HCOOH) Alkalinity (as NH ₃) Carbonyl compounds (as CO) Ethanol (GC) Substances darkened by sulfuric acid Substances reducing permanganate (as O) Sp. Gr.(20/20°C) (Ph.Eur.)	Colorless clear liquid min.99.8%(mass/mass) to pass test 0.791~0.793g/ml 1.327~1.330 max.0.1%(mass/mass) max.5ppm(mass/mass) max.0.002%(mass/mass) max.3ppm(mass/mass) max.0.005%(mass/mass) max.0.05%(mass/mass) to pass test max.5ppm(mass/mass) 0.791~0.793

JIS K 8891:2006*:Additional test performed by Wako but not required by JIS