

SPECIFICATION

Nitric Acid
Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	Colorless clear liquid
Assay (HNO ₃)	65~66%(mass/mass)
Residue after ignition (as sulfate)	max.3ppm(mass/mass)
Chloride (Cl)	max.0.3ppm(mass/mass)
Sulfate (SO ₄)	max.0.5ppm(mass/mass)
Copper (Cu)	max.0.1 ppm(mass/mass)
Lead (Pb)	max.0.1 ppm(mass/mass)
Arsenic (As)	max.0.01ppm(mass/mass)
Iron (Fe)	max.0.2ppm(mass/mass)

JIS K 8541:2015

(*:Additional test performed by Wako but not required by JIS)

The density of this product is approximately 1.40g/mL