

## SPECIFICATION

## Nitric Acid (1.38)

## Super Special Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Chloride (Cl)	max.0.1ppm
Phosphate (PO <sub>4</sub> )	max.0.1ppm
Sulfate (SO <sub>4</sub> )	max.0.1ppm
Sodium (Na)	max.0.1ppm
Potassium (K)	max.0.1ppm
Copper (Cu)	max.0.003ppm
Calcium (Ca)	max.0.01ppm
Zinc (Zn)	max.0.005ppm
Cadmium (Cd)	max.0.1ppb
Mercury (Hg)	max.0.05ppb
Lead (Pb)	max.0.5ppb
Arsenic (As)	max.0.001ppm
Chromium (Cr)	max.0.001ppm
Manganese (Mn)	max.0.001ppm
Iron (Fe)	max.0.01ppm
Cobalt (Co)	max.0.003ppm
Nickel (Ni)	max.0.001ppm
Assay	60.0~62.0%