

Code No.14-0419

## SPECIFICATION

## Nitromethane

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	Colorless ~ slightly yellow, clear liquid
Assay (CH <sub>3</sub> NO <sub>2</sub> )(GC)	min.96.0%(area/area)
Density (20°C)	1.135~1.140g/mL
Refractive index 20°C	1.379~1.384
Water	max.0.2%(mass/mass)
Residue after evaporation	max.0.005%(mass/mass)
Suitability for pesticide analysis	to pass test

JIS K 9523:2019(\*:Additional test performed by Wako but not required by JIS)