

## SPECIFICATION

## 2,2',2''-Nitrilotriethanol

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay [N(CH <sub>2</sub> CH <sub>2</sub> OH) <sub>3</sub> ](GC) Solubility in water Refractive index 20°C Water Residue after ignition (as sulfate) Lead (Pb) Iron (Fe) 2-Aminoethanol (H <sub>2</sub> NCH <sub>2</sub> CH <sub>2</sub> OH) (GC) 2,2'-Iminodiethanol [HN(CH <sub>2</sub> CH <sub>2</sub> OH) <sub>2</sub> ](GC)	Colorless ~ slightly yellow, clear liquid min.98.0%(mass/mass) to pass test 1.482~1.487 max.0.3%(mass/mass) max.0.01%(mass/mass) max.0.001%(mass/mass) max.5ppm(mass/mass) max.0.3%(mass/mass) max.0.3%(mass/mass)

JIS K 8663:2013

(\*:Additional test performed by Wako but not required by JIS)