

Code No.14-0683

## SPECIFICATION

### Nitric Acid, Fuming (1.45)

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	Pale yellow ~ reddish brown, clear liquid
Chloride (Cl)	max.2ppm
Sulfate (SO <sub>4</sub> )	max.5ppm
Heavy metals (as Pb)	max.2ppm
Arsenic (As)	max.0.1ppm
Iron (Fe)	max.1ppm
Assay (as HNO <sub>3</sub> )	77.0~79.0%