

Code No.14-0962

SPECIFICATION

Ninhydrin

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White ~ pale yellow or pale red, crystals ~ crystalline powder to pass test
Solubility in water	
Residue after ignition (as sulfate)	max.0.1%
pH of a 10g/l solution at 25°C	3.5~6.0
Suitability for amino acids determination (E440)	max.0.14
Suitability for amino acids determination (E460)	max.0.04
Suitability for amino acids determination (E570)	max.0.006
Sensitivity	to pass test
Assay	min.98.0%