

Code No.15-0013

SPECIFICATION

1-Octanol

Guaranteed Reagent

| REQUIREMENT | SPECIFICATION |
|---|------------------------|
| * Appearance | Colorless clear liquid |
| Assay [CH ₃ (CH ₂) ₆ CH ₂ OH] (GC) | min.98.0%(mass/mass) |
| Density (20°C) | 0.820~0.833g/ml |
| Refractive index 20°C | 1.428~1.431 |
| Water | max.0.1%(mass/mass) |
| Residue after evaporation | max.0.01%(mass/mass) |
| Acid value | max.0.1 |

JIS K 8213:2014

*:Additional test performed by Wako but not required by JIS