

## SPECIFICATION

## 1-Octanol

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	Colorless clear liquid
Assay [CH <sub>3</sub> (CH <sub>2</sub> ) <sub>6</sub> CH <sub>2</sub> OH] (GC)	min.98.0%(mass/mass)
Density (20°C)	0.820~0.833g/ml
Refractive index 20°C	1.428~1.431
Water	max.0.1%(mass/mass)
Residue after evaporation	max.0.01%(mass/mass)
Acid value	max.0.1

JIS K 8213:2014

\*:Additional test performed by Wako but not required by JIS