

Code No.15-0200

SPECIFICATION

Oxytetracycline Hydrochloride Standard  
for High Performance Liquid Chromatography

| REQUIREMENT                         | SPECIFICATION             |
|-------------------------------------|---------------------------|
| Appearance                          | Yellow, crystals ~ powder |
| Solubility in dil.hydrochloric acid | to pass test              |
| Loss on drying at 105°C             | max.0.5%                  |
| Assay (HPLC)                        | min.99.0%                 |
| Assay (volumetric analysis)         | min.98.0%                 |