

Code No.15-0251

SPECIFICATION

1-Octanol

for Partition Coefficient Analysis

REQUIREMENT	SPECIFICATION
* Appearance	Colorless clear liquid
Assay (GC)	min.99.5%
Density (20°C)	0.824~0.828g/ml
Refractive index 20°C	1.428~1.431
* Absorbance (210nm)	under investigation
* Absorbance (220nm)	under investigation
* Absorbance (230nm)	under investigation
* Absorbance (240nm)	under investigation
* Absorbance (254nm)	under investigation
* Absorbance (279nm)	under investigation
Residue after evaporation	max.0.01%
Acid value	max.0.1

JIS K 8213-1996

*:Additional test performed by Wako but not required by JIS