

Code No.16-0010

## SPECIFICATION

Palmitic Acid  
Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White, small granules or powder
Assay (C <sub>16</sub> H <sub>32</sub> O <sub>2</sub> )(GC)	min.95.0%(mass/mass)
Melting point	60~64°C
Residue after ignition (as sulfate)	max.0.01%(mass/mass)
Acid value	216~221
Iodine value	max.1.5

JIS K 8756:2015

(\*:Additional test performed by Wako but not required by JIS)