

Code No.16-0194

## SPECIFICATION

## Phenyl Salicylate

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White ~ slightly yellow, crystals ~ crystalline powder or mass to pass test
Solubility in ethanol	
Melting point	41~44°C
Chloride (Cl)	max.0.01%
Sulfate (SO <sub>4</sub> )	max.0.01%
Heavy metals (as Pb)	max.0.001%
Phenol and salicylic acid	to pass test
Assay	min.98.0%