

## SPECIFICATION

## Potassium Chromate

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay (K <sub>2</sub> CrO <sub>4</sub> ) Solubility in water pH of a 50g/l solution at 25°C Chloride (Cl) Sulfate (SO <sub>4</sub> ) Sodium (Na) Copper (Cu) Calcium (Ca) Lead (Pb) Iron (Fe)	Yellow, crystals ~ crystalline powder min.99.0%(mass/mass) to pass test 8.5~9.8 max.0.005%(mass/mass) max.0.03%(mass/mass) max.0.05%(mass/mass) max.0.001%(mass/mass) max.0.005%(mass/mass) max.0.001%(mass/mass) max.0.002%(mass/mass)

JIS K 8312:2011

(\*:Additional test performed by Wako but not required by JIS)