

## SPECIFICATION

## Potassium Hexacyanoferrate(II) Trihydrate

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay [K <sub>4</sub> [Fe(CN) <sub>6</sub> ]·3H <sub>2</sub> O] Solubility in water Chloride (Cl) Sulfate (SO <sub>4</sub> ) Sodium (Na)	Pale yellowish brown, crystals ~ crystalline powder min.99.5%(mass/mass) to pass test max.0.003%(mass/mass) max.0.005%(mass/mass) max.0.01%(mass/mass)

JIS K 8802:2007(\*:Additional test performed by Wako but not required by JIS)