

Code No.16-0374

SPECIFICATION

Potassium Hexacyanoferrate(II) Trihydrate

Guaranteed Reagent

| REQUIREMENT | SPECIFICATION |
|-------------------------------------|---|
| * Appearance | Pale yellowish brown, crystals ~ crystalline powder |
| Assay $[K_4[Fe(CN)_6] \cdot 3H_2O]$ | min.99.5%(mass/mass) |
| Solubility in water | to pass test |
| Chloride (Cl) | max.0.003%(mass/mass) |
| Sulfate (SO ₄) | max.0.005%(mass/mass) |
| Sodium (Na) | max.0.01%(mass/mass) |

JIS K 8802:2022(*:Additional test performed by Wako but not required by JIS)