

Code No.16-0406

## SPECIFICATION

### Potassium Nitrite

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White ~ slightly yellow, crystals ~ crystalline powder to pass test
Solubility in water	
Loss on drying at 105°C	max.1.0%
Chloride (Cl)	max.0.005%
Sulfate (SO <sub>4</sub> )	max.0.01%
Sodium (Na)	max.0.5%
Magnesium (Mg)	max.0.002%
Calcium (Ca)	max.0.002%
Lead (Pb)	max.2ppm
Iron (Fe)	max.2ppm
Assay	min.95.0%