

## SPECIFICATION

## Potassium Dihydrogen Phosphate

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White, crystals ~ crystalline powder
Assay (KH <sub>2</sub> PO <sub>4</sub> )	min.99.5%(mass/mass)
Solubility in water	to pass test
Loss on drying at 105°C	max.0.2%(mass/mass)
pH of a 50g/l solution at 25°C	4.2~4.5
Chloride (Cl)	max.5ppm(mass/mass)
Nitrate	to pass test
Sulfate (SO <sub>4</sub> )	max.0.003%(mass/mass)
Heavy metals (as Pb)	max.0.001%(mass/mass)
Sodium (Na)	max.0.02%(mass/mass)
Arsenic (As)	max.1ppm(mass/mass)
Iron (Fe)	max.0.001%(mass/mass)
Ammonium (NH <sub>4</sub> )	max.0.001%(mass/mass)

JIS K 9007:2014 (Amendment 1)

\*:Additional test performed by Wako but not required by JIS