

SPECIFICATION

Tripotassium Phosphate

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	White ~ nearly white, crystals ~ powder
Solubility in water	to pass test
pH of a 50g/L solution at 25°C	min.12.0
Chloride (Cl)	max.0.01%
Sulfate (SO ₄)	max.0.02%
Heavy metals (as Pb)	max.0.002%
Arsenic (As)	max.3ppm
Iron (Fe)	max.0.003%
Assay (after ignition)	min.95.0%