

Code No.16-0435

# SPECIFICATION

## Potassium Diphosphate

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	White powder
Solubility in water	to pass test
pH of a 40g/L solution at 25°C	10.0~11.0
Chloride (Cl)	max.0.005%
Orthophosphate	to pass test
Sulfate (SO <sub>4</sub> )	max.0.01%
Heavy metals (as Pb)	max.0.002%
Iron (Fe)	max.0.002%
Assay	min.95.0%