

SPECIFICATION

Potassium (+)-Tartrate-Water(2/1)

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay [(C ₄ H ₄ K ₂ O ₆) ₂ ·H ₂ O] Solubility in water pH of a 50g/l solution at 25°C Chloride (Cl) Phosphate (PO ₄) Sulfur compounds (as SO ₄) Sodium (Na) Copper (Cu) Calcium (Ca) Zinc (Zn) Lead (Pb) Iron (Fe) Ammonium (NH ₄)	White, crystals ~ crystalline powder min.99.0%(mass/mass) to pass test 7.0~9.0 max.0.001%(mass/mass) max.0.002%(mass/mass) max.0.005%(mass/mass) max.0.5%(mass/mass) max.1ppm(mass/mass) max.0.005%(mass/mass) max.2ppm(mass/mass) max.1ppm(mass/mass) max.0.001%(mass/mass) max.0.002%(mass/mass)

JIS K 8535:2013(*:Additional test performed by Wako but not required by JIS)