

## SPECIFICATION

## Propionic Acid

## Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Solubility in water	to pass test
Density (20°C)	0.991~0.997g/ml
Water	max.0.3%
Chloride (Cl)	max.0.001%
Sulfate (SO <sub>4</sub> )	max.0.005%
Copper (Cu)	max.5ppm
Lead (Pb)	max.5ppm
Iron (Fe)	max.5ppm
Assay (capillary GC)	min.98.0%