

Code No.16-0502

# SPECIFICATION

## Propylene Oxide

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Density (20°C)	0.826~0.831g/mL
Refractive index 20°C	1.364~1.368
Water	max.0.2%
Acidity (as C <sub>2</sub> H <sub>5</sub> COOH)	max.0.005%
Peroxides (as H <sub>2</sub> O <sub>2</sub> )	max.5ppm
Assay (capillary GC)	min.99.0%