

SPECIFICATION

Polyethylene Glycol 400

Wako 1st Grade

| REQUIREMENT | SPECIFICATION |
|-----------------------------------|------------------------|
| Appearance | Colorless clear liquid |
| Viscosity (25°C) | 80~110mPa·s |
| Freezing point | 4~8°C |
| Water | max.1.0% |
| pH of a 50g/l solution at 25°C | 4.0~7.0 |
| Acidity (as CH ₃ COOH) | max.0.1% |
| Average molecular weight | 360~440 |