

Code No.16-0910

SPECIFICATION

Polyethylene Glycol 2,000

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	White ~ nearly white, mass
Viscosity (500g/l, 25°C)	35~60mPa·s
Solubility in water	to pass test
Freezing point	45~55°C
Water	max.1.0%
Residue after ignition (as sulfate)	max.0.2%
pH of a 50g/l solution at 25°C	4.0~7.0
Acidity (as CH ₃ COOH)	max.0.1%
Average molecular weight	1,800~2,200