

Code No.16-0912

## SPECIFICATION

## Polyethylene Glycol 6,000

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	White ~ nearly white, flakes
Viscosity (500g/l, 25°C)	200~400mPa·s
Solubility in water	to pass test
Freezing point	50~65°C
Water	max.1.0%
Residue after ignition (as sulfate)	max.0.2%
pH of a 50g/l solution at 25°C	4.5~7.5
Acidity (as CH <sub>3</sub> COOH)	max.0.1%
Average molecular weight	7,300~9,300