

## SPECIFICATION

Potassium Dihydrogenphosphate for Oxidant Analysis  
for Oxidant Analysis

REQUIREMENT	SPECIFICATION
* Appearance	White crystals
Assay	min.99.0%
Solubility in water	to pass test
pH of a 50g/L solution at 25°C	4.2~4.6
Nitrogen compounds (as N)	max.0.001%
Heavy metals (as Pb)	max.0.001%
Suitability for oxidant determination	to pass test

JIS K 8118:2004

(\*:Additional test performed by Wako but not required by JIS)