

SPECIFICATION

PTH-amino Acids Mixture Standard

for Amino Acid Sequence Analysis

REQUIREMENT	SPECIFICATION
Appearance	Lyophilisate
PTH-L-aspartic Acid	18~22nmol/vial
PTH-L-glutamic Acid	18~22nmol/vial
PTH-L-asparagine	18~22nmol/vial
PTH-L-glutamine	18~22nmol/vial
PTH-DL-serine	18~22nmol/vial
PTH-DL-threonine	18~22nmol/vial
PTH-L-histidine Hydrochloride	18~22nmol/vial
PTH-glycine	18~22nmol/vial
PTH-DL- α -alanine	18~22nmol/vial
PTH-L-tyrosine	18~22nmol/vial
PTH-L-arginine Hydrochloride	18~22nmol/vial
PTH- Δ -threonine	18~22nmol/vial
PTH-DL-methionine	18~22nmol/vial
PTH-DL-valine	18~22nmol/vial
PTH-L-proline	18~22nmol/vial
PTH-DL-tryptophan	18~22nmol/vial
PTH-L-phenylalanine	18~22nmol/vial
PTH-L-lysine	18~22nmol/vial
PTH-DL-isoleucine	18~22nmol/vial
PTH-L-leucine	18~22nmol/vial