

Code No.16-1380

SPECIFICATION

Polyethylene Glycol 2,000,000

Wako 1st Grade

REQUIREMENT	SPECIFICATION	
Appearance	White ~ slightly yellow, crystalline powder ~ powder or granules 25~60mPa·s	
Viscosity (0.5%, 25°C)		
Solubility in water		to pass test
Loss on drying at 105°C		max.5.0%
Residue after ignition (as sulfate)		max.2.0%
pH of a 5g/l solution at 25°C		6.0~8.0