

Code No.16-1537

# SPECIFICATION

## Prometryn Standard for Pesticide Residue Analysis

| REQUIREMENT                              | SPECIFICATION            |
|--|--------------------------|
| Appearance                               | White, crystals ~ powder |
| Solubility in acetone                    | to pass test             |
| Melting point                            | 118~122°C                |
| Loss on drying over silica gel in vacuum | max.0.2%                 |
| Assay (capillary GC)                     | min.98.0%                |
| Assay (qNMR)                             | min.98.0%                |