

Code No.16-1574

SPECIFICATION

Polyoxyethylene(20) Sorbitan Monooleate

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ pale yellow, clear ~ slightly muddy liquid
Viscosity (25°C)	400~480mPa·s
Density (20°C)	1.065~1.095g/ml
Refractive index 20°C	under investigation
Water	max.3.0%
Acid value	max.2
Saponification value	45~55
Iodine value	19~24
Hydroxyl value	under investigation