

Code No.16-1791

SPECIFICATION

Polyvinyl Alcohol

Guaranteed Reagent

| REQUIREMENT | SPECIFICATION |
|---|--|
| * Appearance | White ~ slightly yellow, shot ~ powder |
| Solubility in water | to pass test |
| Loss on drying at 105°C | max.10.0%(mass/mass) |
| Residue after ignition (as sulfate) | max.1.0%(mass/mass) |
| Viscosity (4%, 20°C) | 35~55mPa·s |
| Suitability for absorptiometric determination | to pass test |
| Saponification degree | 78~82mol% |

JIS K 9550 : 2020(*:Additional test performed by Wako but not required by JIS)