

## SPECIFICATION

Pyridine, Dehydrated  
for Organic Synthesis

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Solubility in water	to pass test
Water	max.0.005%
Residue after evaporation	max.0.002%
Chloride (Cl)	max.0.001%
Sulfate (SO <sub>4</sub> )	max.0.001%
Copper	to pass test (Cu:max.about 5ppm)
Ammonium (NH <sub>4</sub> )	max.0.001%
Substances reducing permanganate (as O)	max.0.002%
Assay (capillary GC)	min.99.5%