

Code No.16-1934

# SPECIFICATION

## Potassium Benzoate

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White ~ nearly white, crystals ~ powder or mass
Solubility in water	to pass test
pH of a 50g/l solution at 25°C	7.0~9.0
Chloride (Cl)	max.0.002%
Sulfate (SO <sub>4</sub> )	max.0.01%
Heavy metals (as Pb)	max.0.001%
Iron (Fe)	max.0.001%
Assay	min.99.0%