

Code No.16-2111

SPECIFICATION

Polyoxyethylene(20) Sorbitan Monolaurate for Biochemistry

| REQUIREMENT | SPECIFICATION |
|----------------------|--|
| Appearance | Yellowish brown, clear ~ slightly muddy liquid |
| Viscosity (25°C) | 300~500mPa·s |
| Solubility in water | to pass test |
| Water | max.5.0% |
| Protease | to pass test |
| Acid value | max.3 |
| Saponification value | 40~50 |
| Hydroxyl value | 90~110 |